

The Lowzone Gel2 cushion combines cutting edge technology and over 20 years' experience of living in a wheelchair, to provide a new level of user stability and comfort. This revolutionary cushion encompasses the key factors that affect patient comfort, safety and stability while actively managing thermal regulation and integrating postural and pressure management.

Suitable for medium to high risk categories and for clients with temperature regulating issues. This product should be used in accordance with a therapists recommendation or prescription etc. The Lowzone Gel2 is not a substitute for a comprehensive skin care regime as advised by a clinician.

STORAGE AND USE

The Gel2 cushion can be stored and transported in its' yellow transit bag. Once the cushion is in location it should be stored flat at room temperature (approx. 17°C). Allow the cushion to reach room temperature before use, as the Therma-Gel insert needs to acclimatise to the environment.

Note: If the cushion is not stored horizontally the gel may need massaged to create uniform distribution.

Before positioning the cushion for use, please ensure the seat is clear of objects that may puncture, damage or affect the cushion. Place the cushion onto the chair, flat side down, with the printed area at the back.

CLEANING

Wash the covers only. Remove cover and re-zip before placing in washing machine. Machine wash at 50°C on a short cycle. Hang or tumble dry on a cool setting. Cushion base can be washed down with warm soapy water. DO NOT USE BLEACH OR DRY CLEAN.

WARRANTY AND DISPOSAL

The Lowzone Gel2 cushion has a two year warranty, within which time the Helping Hand Company may replace or repair the cushion. The product warranty does not cover misuse including punctures, burns or tear.

Please dispose of the Lowzone Gel2 cushion at an approved Household Waste Recycling Centre.

IMPORTANT NOTICE

The cushion may affect wheelchair stability, consult wheelchair supplier if you have any concerns

Do not leave the cushion outside overnight or in an area of temperature below 5°C (40°F)

Allow the cushion to warm to room temperature (17°C) before use

The Lowzone Gel2 cushion has a weight limit of 24 stone (152.5kg)